



NUOVO METODO DI RANDOM-LIGATION PER LA GENERAZIONE DI COSTRUTTI CAR OTTIMIZZATI PER LE CELLULE NK

Andrea Sarcinelli, PhD Student

Laboratory of Cell and Gene Therapy for Pediatric Tumors
Department of Pediatric Haematology and Oncology
IRCCS Bambino Gesù Children Hospital
Viale San Paolo 15, Rome, Italy



CONGRESSO
NAZIONALE
AIEOP


ROMA, 22-24 Settembre 2025

CENTRO CONGRESSI
UNIVERSITÀ CATTOLICA
DEL SACRO CUORE

Disclosures of Andrea Sarcinelli

Company name	Research support	Employee	Consultant	Stockholder	Speakers bureau	Advisory board	Other
NONE	NONE	NONE	NONE	NONE	NONE	NONE	NONE


CAR-based immunotherapies



- Ex-vivo gene engineering of patient's cells
- Different generations: 1st, 2nd, 3rd
- Expensive Manufacturing
- Toxicities: CRS, ICANS

CAR-based immunotherapies

CAR T-cells




- Higher toxicities
- Autologous or *Allogenic*
- Slow manufacture
- High costs

FDA-approved

CAR T-cell constructs based on T-cells
costimulatory domains: CD3z, CD28, CD8, 41BB

CAR NK-cells




- Less toxicities: NO GvHD
- Allogenic: healthy donors
- Product “Off-the-Shelf”
- CAR-mediated + Innate NK
- Cheaper and Faster

Clinical trials ongoing

Clinical trials use T-cell-optimized CARs
rather than NK-specific constructs.

CAR-design



CAR-NK clinical trials use T-cell-optimized CARs rather than NK-specific constructs

Integration of T-specific co-stimulatory domains into NK-specific CAR constructs


Integration using a *Random-Ligation* approach

CAR-constructs suitable for T and NK cells




Random Ligation

6 different CAR backbones

Transmembrane	Co-receptors
NKp46	CD3z
NKp30	CD3z
NKp44	DAP12
KIR2DS1	DAP12
KIR2DS2	DAP12
KIR3DS1	DAP12



ScFv: 14G2a
Targeting GD2-antigen






- Directional insertion 5' → 3'
- Reduced self-ligation

29 CAR constructs
were successfully
generated

Generation of a Random Library of CAR constructs

In-vitro characterization on T-cells




In-vitro Antitumor Activity was assessed using IncuCyte S3 live-cell imaging system (Sartorius) on 4 different cell lines: SHSY-5Y, IMR-32, D283, 143b.




In-vitro characterization on T-cells


IMR32 endpoint




IMR32 time-course



SHSY5Y endpoint



SHSY5Y time-course




Neuroblastoma

Medulloblastoma


#1 and #5 best CAR constructs on T-cells
selected for NK testing

Neuroblastoma


Osteosarcoma




D283 endpoint




D283 time-course





143b endpoint



143b time-course



In-vitro characterization on NK-cells




In-vitro Antitumor Activity was assessed using IncuCyte S3 live-cell imaging system (Sartorius) on 4 different cell lines: SHSY-5Y, IMR-32, D283, 143b.




In-vitro characterization on NK-cells

STEP 2: Antitumor assays CAR NK-cells (E:T → 1:2)





Neuroblastoma

Medulloblastoma



Neuroblastoma

Osteosarcoma



● Tumor
■ NT
▲ #1
▼ #5

Conclusions

- Efficient method based on *Sfil* RE able to generate a random library of CAR constructs with a high variability (N=29)
- CAR constructs are functionally active in both T and NK cells (#1 and #5 best candidates)
- Co-stimulatory domains efficiently worked *in-vitro* and they can be exploited for CAR design

FUTURE PROPECTIVES

- To assess CAR activity using *in-vivo* experiments with T and NK cells
- Expand the library of CAR constructs
- Validate more co-stimulatory domains to integrate into the CAR design
- Optimization of *Sfil* cloning method

Acknowledgements

Prof. Franco Locatelli

Head of Department of Haematology/Oncology

Francesca Del Bufalo, Giada Del Baldo, Maria Antonietta De Ioris

Prof. Concetta Quintarelli

Head of Unit of Cell and Gene Therapy for Tumors

Biagio De Angelis

Alessio Ottaviani

Simona Manni

Laura Iaffaldano

Marika Guercio

Pamela De Falco

Maria Luisa D'Amore

Mariasole Aurigemma

Maria Cecilia Quadraccia

Alessio Galimi

Jenny Leopoldina Smith

Matilde Sinibaldi

Stefano Di Cecca

Clelia Geraci

Claudia Pinacchio

Gabriele Canciani

Manuela Ceccarelli

Michele Veneruso

Marilena Margiotta



Gabriele De Marco

Flavia Sadun





In-vitro characterization on NK-cells

STEP 1: Potency assays NT NK-cells



Neuroblastoma

Medulloblastoma



Neuroblastoma

Osteosarcoma

E:T ratio

- 4:1
- 2:1
- 1:1
- 1:2
- 1:4
- 1:8
- 1:16
- 0:1